“A-SPINE” Orthopedic Manual Surgical Instruments

(Non-Sterile)

DESCRIPTION

These instructions for use apply to all models of the “A-SPINE” Orthopedic Manual Surgical Instruments (Non-sterile). The “A-SPINE” Orthopedic Manual Surgical Instruments (Non-Sterile) is all non-powered hand-held devices intended for using with Spinal Fixation System and Polymer Disc Spacer System. It assists the surgeon with implanting proper site of the cervical and thoracic-lumbar vertebra. Please refer to our surgical technical manual for all operating methods. These devices are low-risk items and will not cause harm in clinical use according to the technical manual.

Material

The materials of patient contacting components in “A-SPINE” Orthopedic Manual Surgical Instruments (Non-Sterile) included Type 303 Stainless Steel (ASTM A582), Type 304 Stainless Steel (ASTM A276), Type 316 Stainless Steel (ASTM A276), Type 316L Stainless Steel (ASTM A276), Type 420 Stainless Steel (ASTM A276), Type 465 Stainless Steel (ASTM A564), Type 630 Stainless Steel (ASTM A564), and Ti-6Al-4V (ASTM F136).

The materials of non-patient contacting components in “A-SPINE” Orthopedic Manual Surgical Instruments (Non-Sterile) included PTFE, PP, PPS, PPSU, Silicone and Al. These materials are commonly used for medical instruments.

Warning and Cautions

Please refer to the device label to identify single or multiple use devices and components.

- All reusable instruments must be cleaned, sterilized and sterilized before use.
- Do not use steel wool or abrasives for cleaning.
- Neutral or weakly alkaline cleaners are recommended to avoid the choice of detergents and disinfectants containing the following ingredients: strong organic acids, oxidizing acids, strong alkaline solutions, chlorine and iodine.
- Some device materials may develop changes in mechanical, physical or chemical characteristics under conditions of repeated use, cleaning and re-sterilization that may compromise the integrity of the design and/or material leading to diminished safety, performance and/or compliance with relevant specifications.
- The sterilization parameters are only valid for devices that are adequately cleaned; in addition, instrument tray cannot be stacked in the mechanical washer.
- The Immediate-Use (Flash) Sterilization is not recommended, and should only be used in an emergency for individual instrument.
- The following parameters are valid only when the processes are met with recommendations.
- A-SPINE recommends selecting a mechanical method (washer-disinfector) to clean and disinfect the instrument to ensure its effectiveness and repeatability. Manual methods can be used with an ultrasonic bath only if there is no mechanical method.

Reprocessing

Point of use

- After use (within 2 hours maximum), blood and/or debris is removed from the device to prevent it from drying to the surface.
- It is highly recommended to rinse instruments with running water or transfer instruments to a bath with an aldehyde-free disinfectant solution.

Containment and Transportation

- Soiled devices and non-contaminated devices should be transported separately to avoid contamination.

Preparation for cleaning

- It is recommended to clean contaminated equipment as soon as possible after use.
- The A-SPINE instruments must be cleaned separately from the instrument case and tray.
- The assembly unit refers to the technical guide or other supplemental information to remove the handle/shaft/disassembly kit for cleaning.

Manual: pre-cleaning

- Do not use metal brushes or steel wool for cleaning.
- Soak and/or rinse the instrument into loose dry soil or debris for subsequent removal.
- To ensure that all surfaces are wetted, use a syringe to flush the lumen of the instrument.
- Use a suitable brush to remove dirt, paying particular attention to rough surfaces, threaded areas and joints of device that may contaminate or shield the cleaning process.
- Use a brush of sufficient length and size to clean the lumens with fresh neutral pH enzyme cleaner at least 5 times.
- Soak the instrument into neutral pH enzyme cleaner or detergent which follows the solution manufacturer's instructions for correct exposure time, temperature, and concentration.
- Flush the instrument with running purified water or deionized, and use a syringe to flush the lumen of the instrument.
- Visually inspect any remaining soil and repeat the above steps as necessary.
- Allow to drain on absorbent paper or transfer immediately to the cleaning step.

Manual: Ultrasonic

- Prepare an ultrasonic bath with a fresh cleaning solution at the concentration and temperature specified by the detergent manufacturer.
- Immerse the instruments into bath and activate ultrasonic with frequency 35 - 40 kHz at least 10 minutes.
- Use a suitable brush to remove dirt, paying particular attention to rough surfaces, threaded areas and joints of device that may contaminate or shield the cleaning process.
- Rinse instruments at least 1 minute in running water until all traces of cleaning solution are removed.
- Repeat the above steps until the visual inspection does not have any residual soil.

Mechanical washer

- Soak and/or rinse the instrument to loose dry soil or debris for subsequent removal.
- To ensure that all surfaces are wetted, use a syringe to flush the lumen of the instrument.
- Use a suitable brush to remove soil, paying particular attention to rough surfaces, threaded areas, lumens and joints of device that may contaminate or shield the cleaning process.

- Soak the instrument into neutral pH enzyme cleaner or detergent which follows the solution manufacturer's instructions for correct exposure time, temperature, and concentration.
- Repeat the above steps until the visual inspection does not have any residual soil.
- Place instruments into the washer-disinfector and the lid of the instrument case must be open.
- The instruments should not be overlapped for cleaning efficiency.
- Process instruments reference to below parameters:

Cycle	Minimum Time (minutes)	Temperature	Type of solution
Pre-cleaning	5	20°C	Deionized Water
Main Wash I	10	40°C	0.2% pH-Neutral detergent
Main Wash II	2	20°C	0.2% pH-Neutral detergent
Neutralization/Rinse	1	20°C	Deionized Water
Drying	40	Air-dried at room temperature(23°C), natural convection	N/A

Inspection

Before preparing for sterilization, all instruments should be inspected:

- Cleanness
- Appearance damage, including deformation, corrosion (rust, pitting), deep scratches, flaking, corrosion and cracks.
- Proper function, including movement of joint and sharpness of cutting part.
- Smooth assembly

Note: The end of device life determinate by wear and damage caused by the use. Do not use faulty or defective instruments, it may interference surgery operation.

Packaging

- Considering the ventilation rate, the instruments should not be over filled into the instrument case.
- A-SPINE has verified that standard set of metal lid/case combinations can be reused for steam sterilization.
- The case/tray/lid configurations should be double wrapped in blue/green cloth.
- The responsibility for verifying and evaluating other combinations will be on the end user.

Sterilization

- A-SPINE's devices are only valid by moist heat sterilization.

Note: Immediate-Use (Flash) Sterilization is not recommended, and should only be used in an emergency for individual instrument.

- Sterilization reference to below parameters:

Note: Exposure time means that the chamber maintains sufficient sterilization temperature and pressure during this period, but does not include heating and pressurizing time.

Sterilization type	Gravity displacement
Temperature	121°C
Exposure time (minutes)	30
Package	Rigid container with double wrapped

- Plastic materials can be affected by multiple high temperatures and pressures and cause color changes or deterioration.
- The dry time will be affected by the material and volume of the package. A-SPINE recommends dry time is 30-60 minutes. End users should review and adjust drying time.

Storage

- After sterilization, the device should be stored in a dry, dust-free place. The shelf life depends on the sterile barrier used, the storage method, and the environmental conditions.
- End users should define sterilized shelf life of reusable devices base on sterile barrier and storage conditions.

SERVICING

Please contact customer service center (+886-2-2926-7088) or European Representative (+49-511-39089530), the authorized distributor or dealer. If need more information, the “Intended for Use”, “Surgical Technique” and “A-SPINE Reprocessing Manual” can be downloaded from A-SPINE official website: <http://www.aspine.com.tw/>

European Authorized Representative

mdi Europa GmbH

Langenhagener Str. 71, 30855 Langenhagen Germany

Tel: +49 511 39089530 / Fax: +49 511 39089539

E-mail: info@mdi-europa.com

Web: www.mdi-europa.com

Manufacturer

A-SPINE Asia Co., Ltd.

20F. No. 80, Section 1, Chenggong Road, Yonghe District, New Taipei City 234634, Taiwan

Tel: (886) 2-2926-7088 / Fax: (886) 2-2926-8807

E-mail: service@aspine.com.tw

Web: www.aspine.com.tw

Regarding the information about device safety and performance, please refer to EUDAMED public website:

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