Redmond Cervical Implants

English [EN](#) I-PDS-008-KE Español [ES](#) I-PDS-008-DES 繁體中文 [TW](#) I-PDS-012-DT
 Deutsch [DE](#) I-PDS-008-DDE Português [PT](#) I-PDS-008-DPT
 Français [FR](#) I-PDS-008-DFR Italiano [IT](#) I-PDS-008-DIT

December 2023 C0101-Z-z14 V9 7

EN Redmond Cervical Implants

IMPORTANT NOTICE
 The users of the Redmond Cervical Implants acknowledge that they have read and agreed to the statement in this instruction, ensure that you are familiar with appropriate surgical technique.

DESCRIPTION
 Redmond Cervical Implants as a stand-alone device use screw to secure the cage, The design combines the functionality of a cervical spacer and the benefits of a cervical plate and utilizes the advantages of titanium and PEEK in one implant. The instruments are designed for implantations of the Redmond Cervical Implants. The method of implanting is described in the Redmond Cervical Implants Surgical Technique.

MATERIAL
 All components of the Redmond Cervical Implants are made of PEEK (ASTM F2020) and Titanium alloy (Ti6Al4V, ASTM F136/ISO 5832-3), which are biocompatible, corrosion resistant, non-toxic under biological conditions and interferes as little as possible with imaging processes such as X-ray imaging and computer tomography.

INTENDED USE
 The Redmond Cervical Implants is intended for anterior cervical discectomy and fusion (ACDF) in cervical spine segments (C2-T1). The implant provides stability to restore height and lordosis, providing fusion for mature bone patients.

INDICATION

- Degenerative disk disease
- Spondylosis and retrolisthesis grade 1
- Spinal stenosis
- Revision operation
- Pseudarthrosis

CONTRAINDICATIONS

- Patients with fever or leukocytosis
- Patients with infections associated with the spine (e.g., spondylodiscitis)
- Patients with a history of material allergy or who tend to react to foreign bodies
- Patients with inadequate bone quality or quantity (e.g., severe osteoporosis, osteopenia, osteomyelitis)
- Spinal fracture
- Spinal tumor

WARNING AND PRECAUTIONS

- The surgeon must be thoroughly knowledgeable of the mechanical and limitations of the Redmond Cervical Implants. Adequate training, careful reading of the technical manual and experience in the surgical technique is advised before using the Redmond Cervical Implants.
- Appropriate size selection, material, shape, placement, and fixation play a critical role in implant success. The implant being properly selected to minimize the potential risk of implant failure to the least, thus, close approximations of the size, shape, material and strength of implants to the situation of the individual patient renders maximum result to the implant success. In addition, appropriate placement of implant can reduce the risk to damage soft tissue.
- Patients should be informed of the following reason to cause the possible adverse effects: obesity, coronary heart disease, pregnant woman, above Grade 2 Spondylosis without reduction operation, whole body or nerve-disease, serious osteoporosis, cartilage, steroids treatment necessary or abuse medication.
- Internal factors can improve the stability of the spine, but it must be noted that there is enough space for screw implantation in the vertebral body.
- It is not allowed to be removed after the implant contacts physical tissue or body fluid. The reused product will be caused the serious cross-infection.
- X-ray guidance or other image guiding system is recommended to facilitate precise positioning, reduction in order to avoid or minimize the risk of complications, and check implant position.
- Patients receiving implant should be instructed in full detail of the limitations of the implants, including but not limited to the impact of excessive loading through patients' weight or activity, and should be taught to govern their activities accordingly. Patients should clearly understand that the implanted device is not as strong as normal health bone and will deform, loose, or break when excessive loads are applied on. A patient who is active, debilitated or demented and cannot use weight supporting device properly may be at risk particularly during postoperative rehabilitation period.
- If the patient experiences abnormal pain or other clinical symptoms after the operation, please return for examination.

MRI INFORMATION

The device is classified as "MR conditional" according to ASTM standard F2503. The MR conditional components were tested according to the ASTM standards: F2052, F2182, F2213 and F2119. A patient implanted with the device can undergo an MRI examination under the following conditions:

- Field strengths of 1.5T and 3.0T
- Highest field gradient of 30 T/m (3000 G/cm) or less
- Maximum specific absorption rate (SAR) of 2 W/kg for normal operating mode for a scan time of 30 minutes.
- At 1.5T, the nearest part of the implant shall be at least 20 cm out of the isocenter, and the scan duration shall be at most 15 minutes.

- CAUTION:**
- In the experimental test, heating of a maximum of 2.1°C was measured at 1.5T, 15-minute scan time, SAR 2 W/kg, 20 cm from equiplanar point.
 - In the experimental test, heating of a maximum of 4.9°C was measured at 3.0T, 15-minute scan time, SAR 2 W/kg.

Under these scan conditions, the risks to the patient during an examination are low. To minimize heating, the scan time should be as brief as possible and the SAR should be kept as low as possible. Artifacts: MR imaging in the area of the implants can be impaired by artifacts. In the experimental test, there were artifacts up to 24 mm radially around the implants. The attending physician should conduct a careful risk/benefit assessment.

POSSIBLE SIDE EFFECTS

The following are specific adverse effects which should be understood by the surgeon and explained to the patient. These do not include all adverse effects which can occur with surgery in general, but are important considerations particular to medical internal fixation devices. General surgical risks should be explained to the patient prior to surgery.

- Dural leak
- Nerve damage due to surgical trauma
- Infection
- Pain, discomfort or abnormal sensations due to presence of the device
- Sensitivity or allergic reaction to a foreign body

Bending or fracture of the implant, loosening of the implant

- Delayed union or nonunion
- Decrease in bone density due to stress shielding
- Bursitis

CLEANING & STERILIZATION

- Sterile, locking caps and instruments should be cleaned by ultrasonic cleaner with fresh cleaning solution for 10 minutes. Neutral or weakly alkaline cleaners are recommended to avoid the choice of detergents and disinfectants containing the following ingredients: strong organic acids, oxidizing acids, strong alkaline solutions, chlorine and iodine.
- Avoid the use of corrosive products and/or instruments including abrasive sponges and metal brushes.
- All instruments are in operative condition prior to sterilization.
- The locking cap, screws and instruments are recommended to be sterilized with the A-SPINE designed trays and/or cases.
- Instruction for sterilization: the implants and instruments should be sterilized by steam autoclave following the instructions of the sterilizer manufacturer according to the type of sterilizer used and the method in accordance with the internal hospital guidelines to achieve the degree of sterility of 10⁻⁶. The suggested parameters are as follow:
- Steam Wrapped Gravity Cycle at 121 °C/ 250 °F for 30 minutes, dry time for 30 minutes.
- Steam Wrapped Dynamic-Air-Removal (Pre-vacuum) Cycle at 132 °C/ 270 °F for 4 minutes, dry time for 30 minutes.

PACKAGE, LABELING & STORAGE

- The cage is supplied STERILE, and is sterilized by gamma radiation at least 25 kGy dose. Prior to use, check the product expiration date and verify the integrity of the sterile packaging. Do not use if the package is damaged.
- The locking cap, screws and instruments are supplied NON-STERILE.
- The implants are delivered in individual packages. The implant, package, and insert must be intact at the time of receipt. All legal information required for this type of implant is given in the labeling and the insert of each package.
- Be careful in handling and storage of implants. Cutting, aggressive bending, or scratching of implant surface can significantly reduce the strength and fatigue resistance of the implants. This may cause cracks or non-visible internal stresses that lead to fracture of the implants.
- Implants and instruments should be stored away from the corrosive environments such as salt air, moisture.
- Routine inspection and trial assembly of the instruments and implants are recommended to verify if damage occurred in storage or prior surgery.

SERVICE

When the following conditions occur, please contact A-SPINE, the authorized distributor or dealer. The information should include product description, reference number, lot number, your full contact information, and the type of incident, an accurate description of the incident and consequences, and other related technical information to assist future investigation. Please indicate if a written report is required from A-SPINE, the authorized distributor or dealer.

- In case of claim of dissatisfaction of the identification, reliability, safety, efficacy, performance, or any defect
- Upon receipt of implant, implant defect or damage to the package or insert
- Any serious side effect on patient health or safety, or any life threatening issue or death

CE 2185

EC REP

European Authorized Representative
 mit Europa GmbH
 Langenhager Str. 71, 30855 Langenhagen Germany
 Tel: +49 511 39089530 / Fax: +49 511 39089539
 E-mail: info@nd-europe.com
 Web: www.nd-europe.com

Manufacturer

A-SPINE

A-SPINE Asia Co., Ltd.
 20F., No. 80, Section 1, Chengong Road, Yonghe District,
 New Taipei City 234634, Taiwan
 Tel: (886) 2-2926-7088 / Fax: (886) 2-2926-8807
 E-mail: service@aspinetw.com.tw
 Web: www.aspinetw.com.tw
 Regarding the information about device safety and performance, please refer to EUADAMED public website:
<https://ec.europa.eu/health/euadamed/screen/home>

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