A-Mesh Spinal Titanium Mesh

English 繁體中文 简体中文
I-VBR-001-GE I-VBR-001-DT I-VBR-001-ES

April 2024 C0101-Z-z15 V7 5

EN A-Mesh Spinal Titanium Mesh

IMPORTANT NOTICE
The users of the A-Mesh Spinal Titanium Mesh acknowledge that they have read and agreed to the statement in this insert, which is considered contractual.

DESCRIPTION
The A-Mesh Spinal Titanium Mesh is designed to restore biomechanical integrity throughout the thoracic and lumbar spine following vertebrectomy or corpectomy for patients with spine tumors or fractures. The system provides anterior, middle and posterior vertebral column support both immediately after surgery and for prolonged periods in the absence of bone fusion.

MATERIAL
The A-Mesh Spinal Titanium Mesh is made of Titanium alloy (Ti6Al4V, ASTM F136/ISO 5832-3).

INDICATIONS
The A-Mesh Spinal Titanium Mesh is indicated for use in the thoracolumbar spine (i.e., T1 to L5) to replace a diseased vertebral body resected or excised for the treatment of tumors, to achieve anterior decompression of the spinal cord and neural tissues, and to restore the height of a collapsed vertebral body. The A-Mesh Spinal Titanium Mesh is also indicated for use in vertebral fracture or unstable.

CONTRAINDICATIONS
Contraindications include, but not limited to patients with:

- Severe osteoporosis
- Vertebral fracture or disease with any active or suspected latent infection

1

- Conditions that may place excessive stresses on bone and implants, such as severe obesity or degenerative diseases

POSSIBLE ADVERSE EFFECTS

- Spinal surgery complications: genito-urinary system functional syndrome, gastro-intestinal system functional syndrome, cardiovascular circulatory functional syndrome (including thrombosis), bronchopulmonary syndrome (including pulmonary embolism), bursitis, hematorrhea, myocardial infarction, infection, paralysis or death
- Operation resulting in pachymeninx damage, neuropathogenesis, nerve damage, disease position ossification, vessel and nerve obstacle
- Pain, discomfort, or abnormal sensations due to the presence of the device
- Due to improper implantation, inflammation, early loading and trauma result in loosening, bending, fracture or dislocation
- Delayed or non union of the grafted mass and pseudoarthrosis
- Allergic reactions to implant material
- Decrease in bone density due to bone absorption or stress shielding
- Degenerative changes or instability in segments adjacent to fused vertebral levels

WARNING

- The implants must be implanted only by surgeons having undergone the necessary training in spine surgery the implantation of which must be in accordance with the medical and surgical indications, the potential risks and limitations related to this type of surgery, the contraindications, precautions, and complications.
- Prior to clinical use of these implants, the operating surgeons should thoroughly understand all aspects of the surgical procedure and limitations of the spinal fixation system, The A-Mesh Spinal Titanium Mesh is limited for use by surgeons familiar with the pre- and post-operative treatment, surgical technique, cautions and potential risks related with such spinal surgery. The sufficiency of knowledge of surgical technique, adequate reduction,

2

selection and placement of implants, pre and post-operative patient management are considered essential factors to the delivery of a successful surgical treatment.

- Appropriate size selection, material, shape, placement, and fixation play critical role in implant service life that the decision of selecting proper implant is very important. Implant being properly selected minimize potential risk of implant failure to the least, thus, close approximation of the size, shape, material, strength of implants to the skeleton of the individual patient renders maximum result to the implant service life. Surgical implants must never be reused. An explanted implant should never be re-implanted.
- Great care must be taken to protect the implants from mars, nicks and notches generated by the contact with other metal or abrasive objects. Such injuries produce defects in surface finish of the implants as well as internal stresses which in terms, may become the focal point of implant failure.
- Implant must not be reused. Discard all implants immediately after removal of implants from human body, or mishandled intra-operatively.
- The patient should be warned that non-compliance with postoperative instructions could lead to failure of the device and the possible need thereafter for additional surgery to remove the device.

INSTRUMENTS CLEANING

- All implants and instruments should be cleaned by ultrasonic cleaner with distilled water for 5 minutes. CAUTION: The use of sodium hydroxide (NaOH) is prohibited.
- Avoid the use of corrosive products and/or instruments including abrasive sponges and metal brushes.
- Verify that all instruments are in operative condition prior to sterilization.

INSTRUMENTS FOR STERILIZATION
All implants and instruments should be sterilized by steam autoclave following the instructions of the sterilizer manufacturer according to the type of sterilizer used and the method in accordance with the internal hospital guidelines to achieve the degree of sterility of 10⁴. The suggested parameters are as follow:

- Steam Wrapped Gravity Cycle at 121 °C/ 250 °F for 30 minutes.

PACKAGE, LABELING & STORAGE

- The implants of the A-Mesh Spinal Titanium Mesh are supplied non-sterile.
- The implants and instruments should be sterilized via steam under 121 °C/ 250 °F, 20 PSI for 30 min before using.
- The implants are delivered in individual packages. The implant, package, and insert must be intact at the time of receipt. All legal information required for this type of implant is given in the labeling and the insert of each package.
- The implants may be supplied as a complete set in specially designed trays or boxes, which can be sterilized directly.
- Use care in handling and storage of implants. Cutting, aggressive bending, or scratching of implant surface can significantly reduce the strength & fatigue resistance of the implants. This may cause cracks or non-visible internal stresses that lead to fracture of the implants.
- Implants and instruments should be stored away from the corrosive environments such as salt air, moisture.
- Routine inspection and trial assembly of the instruments and implants are recommended to verify if damage occurred in storage or prior surgery.

4

AFTER SERVICE

- In case of claim of dissatisfaction of the identification, reliability, safety, efficacy, or performance of the A-Mesh Spinal Titanium Mesh, or the services provided by A-SPINE Asia Co., Ltd., please contact A-SPINE, the authorized distributor or dealer.
- Any defect or suspected defect in anyone of the implants should be reported to A-SPINE, the authorized distributor or dealer.
- Upon receipt of implant defect or damage to the package or insert should be reported to A-SPINE, the authorized distributor or dealer for replacement.
- Any serious side effect on patient health or safety, or any life threatening issue or death caused by or associated with the use of the A-Mesh Spinal Titanium Mesh due to component malfunction, damage, or inappropriate direction, should be reported immediately to A-SPINE, the authorized distributor or dealer by fax, e-mail or phone.
- Any claim or report listed above should contain product description, reference number, lot number, your full contact information, and the type of incident, an accurate description of the incident and consequences, and other related technical information to assist future investigation. Please indicate if a written report is required from A-SPINE, the authorized distributor or dealer.

Additional Information regarding the specific implants or instruments may be obtained from:

A-SPINE Asia Co., Ltd.
20F., No. 80, Section 1, Chenggong Road, Yonghe District, New Taipei City 234634, Taiwan
Tel: (886) 2-2926-7088 / Fax: (886) 2-2926-8807
E-mail: service@aspine.com.tw
Web: www.aspine.com.tw

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