SmartLoc Omega Spinal Fixation System Instructions for Use

English 繁體中文
I-SFS-012-GE I-SFS-010-FT

July 2024 C0101-Z-z16 V9

EN SmartLoc Omega Spinal Fixation System Instructions for Use

IMPORTANT NOTICE
The users of the SmartLoc Omega Spinal Fixation System acknowledge that they have read and agreed to the statement in this insert, which is considered contractual.

DESCRIPTION
The SmartLoc Omega Spinal Fixation System is a spinal fixation system composed of pedicle screws, smooth rods, transverse link assemblies, connectors for spinal surgery. Various sizes and styles of the implants are available for optimal adaptation of the pathology of the individual patient. Specific instruments have been designed for use with these implants. For the purpose, function and method of these instruments, please refer to the detail instruction of the SmartLoc Omega Spinal Fixation System Surgical Technique.

MATERIAL
The SmartLoc Omega Spinal Fixation System has different rod materials. These materials include Titanium alloy (Ti6Al4V, ASTM F136/ISO 5832-3), Cobalt-Chromium alloy (ASTM F1537/ISO 5832-12), and PEEK (ASTM F2026) with Tantalum (ASTM F560/ISO 13782). All components aside from the rods are made of Titanium alloy (Ti6Al4V, ASTM F136/ISO 5832-3).

INTENDED USE
The SmartLoc Omega Spinal Fixation System is a pedicle screw fixation intended to use (in the thoracic, lumbar, and sacral spine segments) in skeletally mature patients as an adjunct to fusion for immobilization and stabilization, who present acute or chronic instabilities or deformities.

INDICATIONS

- Degenerative disc disease
- Spondylolisthesis
- Trauma (i.e., fracture or dislocation)
- Spinal stenosis
- Deformities or curvatures (i.e., scoliosis, kyphosis, and/or lordosis)
- Tumor
- Pseudoarthrosis
- Failed previous fusion

In osteoporotic patients, the utilization of cannulated screws, cannulated long-arm screws, perforated screws and/or perforated long-arm screws, in combination with bone cement fixation, is appropriate for the described indications.

CONTRAINDICATIONS

- Any active or suspected latent infection of the spine
- Bone stock abnormalities, or deficiency which can not provide adequate support and/or fixation to the implants
- Pathological obesity
- Open wounds
- Severe osteophyte
- Pregnancy
- Patients with a history of material allergy or who tend to react to foreign bodies
- In osteoporotic patients, screws were utilized without the application of cement.

WARNING AND PRECAUTIONS

- The implants must be implanted only by surgeons having undergone the necessary training in spine surgery the implantation of which must be in accordance with the medical and surgical indications, the potential risks and limitations related to this type of surgery, the contraindications, precautions, and complications.
- Prior to clinical use of these implants, the operating surgeons should thoroughly understand all aspects of the surgical procedure and limitations of the spinal fixation system. The SmartLoc Omega Spinal Fixation System is limited for use by surgeons familiar with the pre- and post-operative treatment, surgical technique, cautions and potential risks related with such spinal surgery. The sufficiency of knowledge of surgical technique, adequate reduction, selection and placement of implants, pre- and post-operative patient management are considered essential factors to the delivery of a successful surgical treatment.
- Appropriate size selection, material, shape, placement, and fixation play critical role in implant service life that the decision of selecting proper implant is very important. Implant being properly selected minimize potential risk of implant failure to the least, thus, close approximation of the size, shape, material, strength of implants to the skeleton of the individual patient renders maximum result to the implant service life.
- The metallic implant does not bear equal strength to resist stress as healthy human bone does. Until maturation of the grafted fusion mass confirmed, implant must not be subject to full weight bearing without additional supports, or implant failure may occur.
- The implant service life may be shortened in case of delay union or nonunion due to metal fatigue. Any scratch, or wear during surgical implantation may become the factor that lead to acceleration of metal fatigue.
- Great care must be taken to protect the implants from mars, nicks and notches generated by the contact with other metal or abrasive objects. Such injuries produce defects in surface finish of the implants as well as internal stresses which in terms, may become the focal point of implant failure.
- Contact of dissimilar metals may accelerate metal corrosion process that occurs on all implanted metals and alloys. The occurrence of corrosion may accelerate fatigue fracture of implants and the increase of metal compound release into human body. Thus, the material of internal fixation such as screws, hooks, connectors, and rods, which may come into contact with another metallic object, should be made of same or compatible

- metals. It is suggested that the SmartLoc Omega Spinal Fixation System not to be used with another system from different sources, manufacturers, or material. If such mixing occurs, A-SPINE claims no responsibility.
- Implant must not be reused. Discard all implants immediately after removal of implants from human body, or mishandled intra-operatively. The reused product will be caused the serious cross-infection.
 - X-ray guidance or other image guided system is recommended to facilitate precise positioning, and reduction in order to avoid or minimize the risk of complications.
 - The SmartLoc Omega Spinal Fixation System consists of implants such as screws, rods, as well as a dedicated set of instruments. Unless otherwise indicated, other instruments should not be used with this system for they may not be compatible the result of which will not be guaranteed.
 - Smoking patients have been shown to exhibit an increased incidence of non-union. Such patients should be warned of this fact as well as its potential consequence.
 - In case that the involvement of patient's occupation or activity which applies excessive load to the implant, failure of implant may occur
 - In some rare cases, progression of degenerative disease may be so advanced that the expected service life of such device may be substantially decreased. In such case, orthopaedic devices may be considered only as a delaying technique or to provide temporary relief.
 - Patients receiving implant should be instructed in full detail of the limitations of the implants, including but not limited to the impact of excessive loading through patients' weight or activity, and should be taught to govern their activities accordingly. Patients should clearly understand that the implanted device is not as strong as normal healthy bone and will bend, deform, loose, or break when excessive loads are applied on. A patient who is active, debilitated or demented and cannot use weight supporting device properly may be at risk particularly during postoperative rehabilitation period.
 - Removing an unlocked spinal screw may require the use of special instrument to disrupt the interface between implant surface and bone tissue. Prior surgery practice of the use of this technique may be required.
 - Prior to making decision of implant removal, surgeon should take into consideration the risk of additional surgery to the patient and the difficulty of removal.

WARNING

The safety and effectiveness of pedicle screw fixation system have been confirmed only under the conditions of deformity or significant instability of the spine requiring fusion with instrumentation. These conditions are: significant instability or deformity of the thoracic, lumbar and sacral spine caused by degenerative spondylolisthesis with distinctive evidence of neurological impairment, scoliosis, fracture, dislocation, spinal tumor, kyphotic deformity, and failed previous fusion. The safety and effectiveness of these devices for other condition remain unknown.

POSSIBLE SIDE EFFECTS

Patients should be advised prior to implantation surgery the possible side effects of such surgeries including:

- Bleeding or hematoma.
- Infection.
- Neurological complications, paralysis, soft tissue lesion, breakage, deformation, or migration of the implants.
- Pain or/and abnormal sensation due to implants.
- Bending, deformation, migration, fracture, or loosening of implants.
- Foreign body allergic reaction to implants.
- Presence of metallic micro particle around the implants.
- Non-union, delayed union, mal-union.
- Bursitis, skin ulcers and or other complications.
- Failure to thoroughly remove the intervertebral disc or properly reposition the sternocleidomastoid muscle can lead to nerve root compression or injury.

Possible side effects including but not limited to the above, may necessitate further surgical treatment.

CLEANING & STERILIZATION

- All implants (except PEEK Rods) and instruments should be cleaned by ultrasonic cleaner with stillled water for 5 minutes.
CAUTION: The use of sodium hydroxide (NaOH) is prohibited.
- Avoid the use of corrosive products and/or instruments including abrasive sponges and metal brushes.
- Verify that all instruments are in operative condition prior to sterilization.
- Instruction for sterilization: The implants and instruments should be sterilized by steam autoclave following the instructions of the sterilizer manufacturer according to the type of sterilizer used and the method in accordance with the internal hospital guidelines to achieve the degree of sterility of 10-6. The suggested parameters are as follow:
 - Steam Wrapped Gravity Cycle at 121 ℃/ 250 ℉ for 30 minutes, dry time for 30 minutes.
 - Steam Wrapped Dynamic-Air-Removal (Pre-vacuum) Cycle at 132 ℃/ 270 ℉ for 4 minutes, dry time for 30 minutes.

After processing, please inspect the cleaning, damage and function of implants and instruments prior to use.

PACKAGE, LABELING & STORAGE

- The PEEK Rods had been sterilized by gamma radiation. It should avoid contaminating while operation process. It is necessary to exchange if the packaging has been broken without reason.
- The other implants and instruments are supplied NON-STERILE.
- The implants are delivered in individual packages. The implant, package, and insert must be intact at the time of receipt. All legal information required for this type of implant is given in the labeling and the insert of each package.
- The implants may be supplied as a complete set in specially designed trays or boxes, which can be sterilized directly.

- Use care in handling and storage of implants. Cutting, aggressive bending, or scratching of implant surface can significantly reduce the strength & fatigue resistance of the implants. This may cause cracks or non-visible internal stresses that lead to fracture of the implants.
- Implants and instruments should be stored away from the corrosive environments such as salt air, moisture.
- Routine inspection and trial assembly of the instruments and implants are recommended to verify if damage occurred in storage or prior surgery.

SERVICE

When the following conditions occur, please contact A-SPINE, the authorized distributor or dealer. The information should include product description, reference number, lot number, your full contact information, and the type of incident, an accurate description of the incident and consequences, and other related technical information to assist future investigation. Please indicate if a written report is required from A-SPINE, the authorized distributor or dealer.

- In case of claim of dissatisfaction of the identification, reliability, safety, efficacy, performance, or any defect of the SmartLoc Omega Spinal Fixation System
- Upon receipt of implant, implant defect or damage to the package or insert
- Any serious side effect on patient health or safety, or any life threatening issue or death

A-SPINE Asia Co., Ltd.
20F., No. 80, Section 1, Chenggong Road, Yonghe District, New Taipei City 234634, Taiwan
Tel: (886) 2-2926-7088 / Fax: (886) 2-2926-8807
E-mail: service@aspine.com.tw
Web: www.aspine.com.tw

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