Polymer Disc Spacer
“A-Spine” Vigor PEEK Cervical Disc Spacer

English 繁體中文 简体中文
I-PDS-003-GE I-PDS-001-DT I-PDS-001-FS

January
2024

C0101-Z-z09
V10

EN

Polymer Disc Spacer
“A-Spine” Vigor PEEK Cervical Disc Spacer

DESCRIPTION

The “A-Spine” Vigor PEEK Cervical Disc Spacer is made of PEEK (ASTM F2026), Tantalum (ASTM F560/ ISO 13782) and Titanium alloy (Ti6Al4V, ASTM F136/ISO 5832-3). The instruments are designed for implantations of the “A-Spine” Vigor PEEK Cervical Disc Spacer. The method of implanting is described in the “A-Spine” Vigor PEEK Cervical Disc Spacer Surgical Technical Manual.

INDICATIONS

Indications for fusion procedures are for example:

- Degenerative disk disease
- Spondylolisthesis and retrolisthesis grade 1
- Degenerative scoliosis
- Revision operation
- Pseudarthrosis

CONTRAINDICATIONS

- Patients with fever or leukocytosis
- Patients with infections associated with the spine (e.g. spondylodiscitis)
- Patients with a history of material allergy or who tend to react to foreign bodies
- Patients whose general medical or psychological condition is unfavorable for- or could be worsened by the procedure; careful consideration is required on the part of the treating physician/surgeon for these patients
- Patients with inadequate bone quality or quantity (e.g. severe osteoporosis, osteopenia, osteomyelitis)
- Pregnancy

PRECAUTIONS

- The surgeon must be thoroughly knowledgeable of the mechanical and limitations of the “A-Spine” Vigor PEEK Cervical Disc Spacer. Adequate training, careful reading of the technical manual and experience in the surgical technique is advised before using the “A-Spine” Vigor PEEK Cervical Disc Spacer.

- Patients should be informed of the following reason to cause the possible adverse effect: Obesity, coronary heart disease, pregnant woman, above Grade 2
- The surgeon must be thoroughly knowledgeable of the mechanical and limitations of the “A-Spine” Vigor PEEK Cervical Disc Spacer. Adequate training, careful reading of the technical manual and experience in the surgical technique is advised before using the “A-Spine” Vigor PEEK Cervical Disc Spacer.
- Patients should be informed of the following reason to cause the possible adverse effects: Obesity, coronary heart disease, pregnant woman, above Grade 2 Spondylolisthesis without reduction operation, whole body or nerve-end disease, serious osteoporosis, cartilage, steroids treatment necessarily or abuse medicament.
- It is better to implant internal fixator increasing the stability of spine.
- It is not allowed to be re-used after the implant contacts physical tissue or body fluid. The reused product will be caused the serious cross-infection.

POSSIBLE ADVERSE EFFECTS

The following things are specific adverse effects which should be understood by the surgeon and explained to the patient. These do not include all adverse effects which can occur with surgery in general, but they are important considerations particular to metallic internal fixation devices. General surgical risks should be explained to the patient prior to surgery.

- Dural leak
- Nerve damage due to surgical trauma
- Infection
- Pain, discomfort or abnormal sensations due to presence of the device
- Sensitivity or allergic reaction to a foreign body
- Bending or fracture of the implant; loosening of the implant
- Delayed union or nonunion
- Decrease in bone density due to stress shielding
- Bursitis

STERILIZATION

The cage had been sterilized by gamma radiation at least 25 kGy dose. It should avoid contaminating while operation process. It is necessary to exchange if the packaging has been broken without reason. The instruments are provided Non-Sterile to hospital. A-SPINE suggests the device should be sterilized via steam sterilization (under 121 °C/ 250 °F, 20 PSIG for 30 mins) to assure SAL (10⁻⁶).

SERVICING

Please contact customer service center (+886-2-2926-7088), the authorized distributor or dealer, meanwhile, please offer patients' information and item needed in order to preparation and delivery.

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