Winloc Anterior Cervical Plate System Instructions for Use

English 繁體中文 简体中文
I-SFS-005-HE I-SFS-016-ET I-SFS-007-FS

May 2024 C0101-Z-z10 V10

Winloc Anterior Cervical Plate System Instructions for Use

IMPORTANT NOTICE

The users of the Winloc Anterior Cervical Plate System acknowledge that they have read and agreed to the statement in this instruction, ensure that you are familiar with appropriate surgical technique.

DESCRIPTION

The Winloc Anterior Cervical Plate System consists of wide range of cervical plates and screws. The bone plate design can prevent screws from backing out and can be combined with expandable bone screws, enhancing stability of fixation and leading to improved therapeutic results. The instruments are designed for implantations of the Winloc Anterior Cervical plate system. The method of implanting is described in the Winloc Anterior Cervical Plate System Surgical Technique.

MATERIAL

All components of the Winloc Anterior Cervical plate system are made of Titanium alloy (Ti6Al4V, ASTM F136/ISO 5832-3), which are biocompatible, corrosion resistant, non-toxic under biological conditions and interferes as little as possible with imaging processes such as X-ray imaging and computer tomography.

INTENDED USE

The Winloc Anterior Cervical Plate System is intended for the anterior surgical stabilization of the cervical spine. The implant enables the immobilization of the fusion mass until consolidation of the implant or bone transplantation.

INDICATIONS

- Degenerative disc
- Disease
- Spondylolisthesis
- Trauma (i.e., fracture or dislocation)
- Spinal stenosis
- Deformities or curvatures(i.e., scoliosis, kyphosis,and/or lordosis)
- Tumor
- Pseudoarthrosis
- Failed previous fusion

CONTRAINDICATIONS

- Patients with acute infection, whether superficial or deep
- Bony abnormalities, unable to ensure fixation of bone screws
- Pathological obesity
- Open wounds
- Patients with a history of material allergy or who tend to react to foreign bodies
- The physician must consider carefully before treating patients who are in a generally unfavorable medical or psychological state and who could be made worse by the procedure
- Pregnancy

WARNING AND PRECAUTIONS

- Appropriate size selection, material, shape, placement, and fixation play a critical role in implant service life that the decision of selecting proper implant is very important. Implant being properly selected minimize potential risk of implant failure to the least, thus, close approximation of the size, shape, material, strength of implants to the skeleton of the individual patient renders maximum result to the implant service life.
- The resistance of metallic internal fixation devices differs from that of healthy bones. Hence, it is difficult to anticipate the point of damage for the implant when subjected to body weight without supplementary support.
- Implants may collapse under excessive load, leading to delayed fusion or non-fusion. Therefore, the internal fixator must share the load until the physician verifies bone healing.
- Delayed or non-union healing may result in shortened lifespan of implants due to metal fatigue. The surgeons should avoid notching or scratching the implants which will contribute to early fatigue of the implants.
- Mixing metals or contact between dissimilar metal, such as titanium screws used with a stainless steel bone plate will accelerate the corrosion process of stainless steel. Therefore, devices that come into contact with each other should be made from similar or compatible metals.
- Implant should not be reused and should be discarded upon removal.
- Utilizing X-Ray during surgery for positioning and realignment is crucial in reducing complications.
- The surgeon must be thoroughly knowledgeable of the mechanical and metalurgical limitations of the Winloc Anterior Cervical Plate System implants. Adequate training, careful reading of the technical manual and experiences in the surgical technique is advised before using the Winloc Anterior Cervical Plate System.

- The patient's activity degree will have a significant impact on the useful life span of any temporary internal fixation devices. The Patient must be informed that any activity increases the risk of loosening, bending or breaking the device, especially before proper bone grafting and normal healing of bone tissue occurs.
- Correct handling of the implant is extremely important. Contouring of metal implants should only be performed with proper bending equipment.

POSSIBLE SIDE EFFECTS

The following are specific adverse effects which should be understood by the surgeon and explained to the patient. These do not include all adverse effects which can occur with surgery in general, but are important considerations particular to metallic internal fixation devices. General surgical risks should be explained to the patient prior to surgery.

- The potential complications of spinal surgery may include: urological and reproductive dysfunction, gastrointestinal disturbances, vascular circulatory system dysfunction (including thrombosis), bronchial and pulmonary obstruction (including embolism), Bursitis, haemorrhage, myocardial infarction, infection, paralysis, or death.
- Injury to the dura mater and arachnoid mater caused by surgery. Peripheral nerve abnormalities, nerve injury, localized ossification, vascular and neural disorders.
- Pain, discomfort or abnormal sensations due to presence of the device.
- Improper installation, inflammation, premature loading, and trauma may result in loosening, bending, fracturing, or detachment of the implant.
- Sensitivity or allergic reaction to a foreign body
- Decrease in bone density due to stress shielding

- Delayed union, nonunion or pseudoarthrosis
- Sensitivity in skin or muscle patients could lead to skin ulceration or other complications
- Inadequate excision of intervertebral discs or failure to reposition muscles such as the sternocleidomastoid can cause compression or injury to the nerve roots
- In the case of a serious adverse event, additional surgery might be required

CLEANING & STERILIZATION

- All implants and instruments should be cleaned by ultrasonic cleaner with fresh cleaning solution for 10 minutes. Neutral or weakly alkaline cleaners are recommended to avoid the choice of detergents and disinfectants containing the following ingredients: strong organic acids, oxidizing acids, strong alkaline solutions, chlorine and iodine.
- Avoid the use of corrosive products and/or instruments including abrasive sponges and metal brushes.
- Verify that all instruments are in operative condition prior to sterilization.
- The implants and instruments are recommended to be sterilized with the A-SPINE designed trays and/or cases.
- Instruction for sterilization: the implants and instruments should be sterilized by steam autclave following the instructions of the sterilizer manufacturer according to the type of sterilizer used and the method in accordance with the internal hospital guidelines to achieve the degree of sterility of 10⁻⁴. The suggested parameters are as follow:
 - Steam Wrapped Gravity Cycle at 121 ℃/ 250 ℉ for 30 minutes, dry time for 30 minutes,
 - Steam Wrapped Dynamic-Air-Removal (Pre-vacuum) Cycle at 132 ℃/ 270 ℉ for 4 minutes, dry time for 30 minutes.
- After processing, please inspect the cleaning, damage and function of implants and instruments prior to use.

PACKAGE, LABELING & STORAGE

- The implants and instruments are supplied NON-STERILE.
- The implants are delivered in individual packages. The implant, package, and insert must be intact at the time of receipt. All legal information required for this type of implant is given in the labeling and the insert of each package.
- The implants may be supplied as a complete set in specially designed trays or boxes, which can be sterilized directly.
- Be careful in handling and storage of implants. Cutting, aggressive bending, or scratching of implant surface can significantly reduce the strength & fatigue resistance of the implants. This may cause cracks or non-visible internal stresses that lead to fracture of the implants.
- Implants and instruments should be stored away from the corrosive environments such as salt air, moisture.
- Routine inspection and trial assembly of the instruments and implants are recommended to verify if damage occurred in storage or prior surgery.

SERVICE

When the following conditions occur, please contact A-SPINE, the authorized distributor or dealer. The information should include product description, reference number, lot number, your full contact information, and the type of incident, an accurate description of the incident and consequences, and other related technical information to assist future investigation. Please indicate if a written report is required from A-SPINE, the authorized distributor or dealer.

- In case of claim of dissatisfaction of the identification, reliability, safety, efficacy, performance, or any defect
- Upon receipt of implant, implant defect or damage to the package or insert
- Any serious side effect on patient health or safety, or any life threatening issue or death

A-SPINE Asia Co., Ltd.
20F., No. 80, Section 1, Chengong Road, Yonghe District, New Taipei City 234634, Taiwan
Tel: (886) 2-2926-7088 / Fax: (886) 2-2926-8807
E-mail: service@aspine.com.tw
Web: www.aspine.com.tw

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