Polymer Disc Spacer X'Plo Polymer Lumbar Cage

English 繁體中文 简体中文
I-PDS-015-FE I-PDS-006-DT I-PDS-003-FS

March 2024 C0101-Z-z06 V9

Polymer Disc Spacer X'Plo Polymer Lumbar Cage

DESCRIPTION
The X'Plo Polymer Lumbar Cage is made of PEEK (ASTM F2026) and Tantalum (ASTM F560/ ISO 13782). The instruments are designed for implantations of the X'Plo Polymer Lumbar Cage. The method of implanting is described in the X'Plo Polymer Lumbar Cage Surgical Technical Manual.

INDICATIONS
The diseases use with autogenous bone graft for spinal interbody fusion operation, including:
• Use for Degenerative Disc Disease (DDD) and Degenerative Lumbar Scoliosis at 1 or 2 levels from L1 to S1, e.g. primary laminectomy for decompression
• Grade 1 spondylolisthesis or retrolisthesis at the involved level(s)
• Revision surgery for failed column operation or post-operation instability
• Stenosis
• Pseudarthrosis at the lumbar
• Posterior or anterior approach for lumbar

CONTRAINDICATIONS
• Patients with fever or leukocytosis
• Patients with infections associated with the spine (e.g. spondylodiscitis)
• Patients with a history of material allergy or who tend to react to foreign bodies
• Patients whose general medical or psychological condition is unfavorable for- or could be worsened by the procedure; careful consideration is required on the part of the treating physician/surgeon for these patients
• Patients with inadequate bone quality or quantity (e.g. severe osteoporosis, osteopenia, osteomyelitis)
• Pregnancy

POSSIBLE ADVERSE EFFECTS
The following are specific adverse effects which should be understood by the surgeon and explained to the patient. These do not include all adverse effects which can occur with surgery in general, but are important considerations particular to metallic internal fixation devices. General surgical risks should be explained to the patient prior to surgery.

- Dural leak
- Nerve damage due to surgical trauma
- Infection
- Pain, discomfort or abnormal sensations due to presence of the device
- Sensitivity or allergic reaction to a foreign body
- Bending or fracture of the implant; loosening of the implant
- Delayed union or nonunion
- Decrease in bone density due to stress shielding
- Bursitis

STERILIZATION
The cage had been sterilized by gamma radiation at least 25 kGy dose. It should avoid contaminating while operation process. It is necessary to exchange if the packaging has been broken without reason. The instruments are provided Non-Sterile to hospital. A-SPINE suggests the device should be sterilized via steam sterilization (under 121 °C/ 250 °F, 20 PSIG for 30 mins) to assure SAL (10⁻⁶).

SERVICING
Please contact customer service center (+886-2-2926-7088), the authorized distributor or dealer, meanwhile, please offer patients' information and item needed in order to preparation and delivery.

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