


Combo® Lumbar Disc Cage (Titanium- PEEK)



Combo® utilizes the advantages of titanium & PEEK in one implant

- ◆ Titanium
 - Direct bony ingrowth increases initial implant stability and fusion rate
 - Easier verification of Implant placement
- ◆ PEEK
 - Vertebra-like elasticity reduces the stress shielding
 - Radiolucency allows for reliable healing site monitoring
- ◆ **Combo® : A unique combination of titanium and PEEK**
 - Encouraging rapid bony ingrowth and increasing fusion rate
 - Optimizing intervertebral stability
 - Excellent visualization on postoperative CT & MRI scan

Specifications

Combo® Lumbar Disc Cage

Item Number	Angulation	H/mm	L/mm
595-07248H	5°	7	24
595-08248H	5°	8	24
595-09248H	5°	9	24
595-10248H	5°	10	24
595-11248H	5°	11	24
595-12248H	5°	12	24
595-13248H	5°	13	24
595-14248H	5°	14	24



Comparison

