


Combo® Lumbar Disc Cage (Titanium - PEEK)



Combo® utilizes the advantages of titanium & PEEK in one implant

- ♦ **Titanium**
 - Direct bony ingrowth increases initial implant stability and fusion rate
 - Easier verification of Implant placement
- ♦ **PEEK**
 - Vertebra-like elasticity reduces the stress shielding
 - Radiolucency allows for reliable healing site monitoring
- ♦ **Combo® : A unique combination of titanium and PEEK**
 - Encouraging rapid bony ingrowth and increasing fusion rate
 - Optimizing intervertebral stability
 - Excellent visualization on postoperative CT & MRI scan

Specifications



Item Number	Angulation	H / mm	L / mm
595-07288H	5°	7	28
595-08288H		8	
595-09288H		9	
595-10288H		10	
595-11288H		11	
595-12288H		12	
595-13288H		13	
595-14288H		14	

Item Number	Angulation	H / mm	L / mm
598-07288H	5°	7	28
598-08288H		8	
598-09288H		9	
598-10288H		10	
598-11288H		11	
598-12288H		12	
598-13288H		13	
598-14288H		14	

Comparison

