


Combo[®] Cervical Disc Cage (Titanium - PEEK)



Combo[®] utilizes the advantages of titanium & PEEK in one implant

- ♦ Titanium
 - Direct bony ingrowth increases initial implant stability and fusion rate
 - Easier verification of Implant placement
- ♦ PEEK
 - Vertebra-like elasticity reduces the stress shielding
 - Radiolucency allows for reliable healing site monitoring
- ♦ **Combo[®] : A unique combination of titanium and PEEK**
 - Encouraging rapid bony ingrowth and increasing fusion rate
 - Optimizing intervertebral stability
 - Excellent visualization on postoperative CT & MRI scan

Specifications

Combo® Cervical Disc Cage

Item Number	H/mm	W/mm	L/mm
165-25128H	5	14	12
165-26128H	6	14	12
165-27128H	7	14	12
165-28128H	8	14	12
165-29128H	9	14	12
165-20128H	10	14	12
165-25148H	5	15	14
165-26148H	6	15	14
165-27148H	7	15	14
165-28148H	8	15	14
165-29148H	9	15	14
165-20148H	10	15	14



Comparison

